

AD-TECH® EPILEPSY/LTM BEHNKE FRIED DEPTH ELECTRODES


- For simultaneous macro and micro contact recordings.
- 8-contact macro recording outer body with 8 or 9-contact micro recording.
- **Platinum** contacts.
- Contact 1 starts 1.0 mm +/- 0.5 mm from body tip.
- Sterile; single-use only.

Catalog #: BF08R-SP05X-000

8-contact macro depth electrode.



Catalog #: BF08R-SP21X-0B0

8-contact macro depth electrode.
1.5mm long contacts, 1 thru 2=3mm, 2 thru 7= 5.5mm spacing includes block with through hole



Catalog #: WB09R-SP00X-014

9-contact micro-wire bundle
(8 recording & 1 reference contact). Ends in "standard" pigtail.



Catalog #	Description ¹	Contact Size ²	TECH-ATTACH Components		Cable
			Cable	Block	
BF08R-SP05X-000	8-contact macro depth	1.57/1.28	L-DCL-8DIN	Included	n/a
WB09R-SP00X-0B6	9-contact micro depth	0.04	L-DCL-9DIN	DC-9X	L-SRL-9DIN
WB09R-SP00X-014	9-contact micro depth	0.04	L-DCL-9DIN	DC-9X	L-SRL-9DIN

¹ Each electrode requires the TECH-ATTACH or wire bundle cable as indicated. A special, slotted connector block that connects the macro contacts to the TECH-ATTACH cable is included with each BF-series electrode. Additional slotted blocks can be purchased, contact your Ad-Tech representative for catalog numbers and prices.

² Overall length of the contact followed by its diameter (mm).

³ Micro wires extend 2 cm past the tip; trim to desired length prior to use.

